

L435AG75



## Aeratore F22x100

### CONFORMIDAD

### ADVANTAGES

- ISO 3864 compliant
- First-aid in case of chemical projection



### CARACTERISTICAS TECNICAS

Modelo

L435AG75

Código Aduanero

73242921

Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Caudal de lavajos (L/min)	Caudal de la ducha de seguridad (L/min)	Peso (kg)
L435AG75			0	0	0.60