

9003ISAI



Mela di doccia in acciaio inossidabile

CONFORMIDAD

ADVANTAGES

- ISO 3864 compliant
- First-aid in case of chemical projection



CARACTERISTICAS TECNICAS

Modelo

9003ISAI

Código Aduanero

73242921

Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Caudal de lavajos (L/min)	Caudal de la ducha de seguridad (L/min)	Peso (kg)
9003ISAI			0	0	0.60