

P992IS



Paratoia fuori gelo per DLI63

CONFORMIDAD

ADVANTAGES

- ISO 3864 compliant
- First-aid in case of chemical projection



CARACTERISTICAS TECNICAS

Modelo	P992IS		Código Aduanero	73242921	
Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Caudal de lavaojos (L/min)	Caudal de la ducha de seguridad (L/min)	Peso (kg)
P992IS			0	0	0.60



138 avenue Henri BARBUSSE 92140 Clamart (France)
Tél. 01 46 45 80 00 - Fax : 01 46 45 85 98
www.trionyx.fr - [trionyx\[at\]trionyx.fr](mailto:trionyx[at]trionyx.fr)