

B235-900



traghetto di ritenzione per armadio 2 porte - L 900 - P 615 mm

CONFORMIDAD

ADVANTAGES

- 30-minute fire resistant and made in France
- Storage of flammables, corrosives and toxics



CARACTERISTICAS TECNICAS

Modelo
Dimensiones de envío alt x anc x pro (mm)

B235-900
1000 x 500 x 500

Código Aduanero

94032080

Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Capacidad de almacenamiento (L)	Cubetas de retención	Estanterías	Carga admisible	Peso (kg)
B235-900	100 x 765 x 410		0		0	0	10.00