

C235



Carabottino per traghettino di ritenzione per armadio 1 porta - L 1135 - P 615 mm

CONFORMIDAD

ADVANTAGES

- 30-minute fire resistant and made in France
- Storage of flammables
- Storage of flammables, corrosives and toxics



CARACTERISTICAS TECNICAS

Modelo
Dimensiones de envío alt x anc x pro (mm)

C235
1000 x 500 x 500

Código Aduanero

94032080

Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Capacidad de almacenamiento (L)	Cubetas de retención	Estanterías	Carga admisible	Peso (kg)
C235	90 x 980 x 400		0		0	0	8.00