

C148



## Carabottino per tragheto di ritenzione per armadio 1 porta - L 635 - P 615 mm

### CONFORMIDAD

### ADVANTAGES

- 30-minute fire resistant and made in France
- Storage of flammables



### CARACTERISTICAS TECNICAS

Modelo  
Dimensiones de envío alt x anc x pro (mm)

**C148**  
**700 x 500 x 500**

Código Aduanero

**94032080**

Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Capacidad de almacenamiento (L)	Cubetas de retención	Estanterías	Carga admisible	Peso (kg)
C148	90 x 460 x 390		0		0	0	4.00