

Etagere/bac di ritenzione supplementare per gli scompartimenti acidi, basi, tossici

CONFORMIDAD

ADVANTAGES

- Storage of corrosives
- Retention equipment
- Built in polypropylene



CARACTERISTICAS TECNICAS

Modelo		E8A2		Código Aduanero		39269098	
Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Capacidad de almacenamiento (L)	Cubetas de retención	Estanterías	Carga admisible	Peso (kg)
E8A2	85 x 455 x 480		0		0	20	5.00