

B235



## traghetto di ritenzione per armadio 2 porte - L 1135 - P 615 mm

### CONFORMIDAD

### ADVANTAGES

- 30-minute fire resistant and made in France
- Storage of flammables, corrosives and toxics



### CARACTERISTICAS TECNICAS

Modelo  
Dimensiones de envío alt x anc x pro  
(mm)

**B235**  
**1000 x 500 x 500**

Código Aduanero

**94032080**

Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Capacidad de almacenamiento (L)	Cubetas de retención	Estanterías	Carga admisible	Peso (kg)
B235	100 x 1000 x 410		0		0	0	10.00