

BR250



PE



## Tavolozza di ritenzione coperta per 2 fusti di 200 L

### CONFORMIDAD

### ADVANTAGES

- Retention equipment
- Built in polyethylene
- Storage of flammables, corrosives and toxics

### SINO TAMBIÉN...

- Available in different sizes and capacities



### CARACTERISTICAS TECNICAS

Modelo	BR250			Código Aduanero	94032090		
Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Capacidad de almacenamiento (L)	Cubetas de retención	Estanterías	Carga admisible	Peso (kg)
BR250	2110 x 1565 x 995		400		0	0	101.00