

BR1000



PE



## Tavolozza di ritenzione per tino 1000 L

### CONFORMIDAD

### ADVANTAGES

- Retention equipment
- Built in polyethylene
- Storage of flammables, corrosives and toxics

### SINO TAMBIÉN...

- Available in different sizes and capacities



### CARACTERISTICAS TECNICAS

| Modelo | BR1000                                      |   |                                 | Código Aduanero      | 94032090    |                 |           |
|--------|---|---|---------------------------------|----------------------|-------------|-----------------|-----------|
| Modelo | Dimensiones exteriores alt x anc x pro (mm) | Dimensiones interiores alt x anc x pro (mm) | Capacidad de almacenamiento (L) | Cubetas de retención | Estanterías | Carga admisible | Peso (kg) |
| BR1000 | 1000 x 1450 x 1450                          |   | 1000                            |                      | 0           | 0               | 105.00    |