

BR1000



PE



## Tavolozza di ritenzione per tino 1000 L


### CONFORMIDAD

### ADVANTAGES

- Retention equipment
- Built in polyethylene
- Storage of flammables, corrosives and toxics

### SINO TAMBIÉN...

- Available in different sizes and capacities



### CARACTERISTICAS TECNICAS

Modelo	BR1000			Código Aduanero	94032090		
Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Capacidad de almacenamiento (L)	Cubetas de retención	Estanterías	Carga admisible	Peso (kg)
BR1000	1000 x 1450 x 1450		1000		0	0	105.00