

28234

NFPA  
OSHA

## Tavolozza di ritenzione per 2 fusti di 200 litri

### CONFORMIDAD

EPA - SPCC - U - O - UL

### ADVANTAGES

- Retention equipment
- NFPA and OSHA compliant
- Storage of flammables, corrosives and toxics

### SINO TAMBIÉN...

- EPA, SPCC, UL Environment, NFPA and OSHA compliant
- Conforme con lo artículo 2 Febrero 1998



### CARACTERISTICAS TECNICAS

Modelo  
Conformidad**28234**  
**EPA - SPCC - U - O - UL**

Código Aduanero

**94032090**

Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Capacidad de almacenamiento (L)	Cubetas de retención	Estanterías	Carga admisible	Peso (kg)
28234	457 x 1245 x 635		400		0	0	20.00