

28703

NFPA
OSHA

Tavolozza di ritenzione per 3 fusti di 200 litri

CONFORMIDAD

EPA - SPCC - U - O - UL

ADVANTAGES

- Retention equipment
- NFPA and OSHA compliant
- Storage of flammables, corrosives and toxics

SINO TAMBIÉN...

- EPA, SPCC, UL Environment, NFPA and OSHA compliant
- Conforme con lo artículo 2 Febrero 1998



CARACTERISTICAS TECNICAS

Modelo
Conformidad**28703**
EPA - SPCC - U - O - UL

Código Aduanero

94032090

Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Capacidad de almacenamiento (L)	Cubetas de retención	Estanterías	Carga admisible	Peso (kg)
28703	295 x 1854 x 635		600		0	0	26.00