

28704

NFPA  
OSHA

## Tavolozza di ritenzione per 4 fusti di 200 litri

### CONFORMIDAD

EPA - SPCC - U - O - UL

### ADVANTAGES

- Retention equipment
- NFPA and OSHA compliant
- Storage of flammables, corrosives and toxics

### SINO TAMBIÉN...

- EPA, SPCC, UL Environment, NFPA and OSHA compliant
- Conforme con lo artículo 2 Febrero 1998



### CARACTERISTICAS TECNICAS

Modelo  
Conformidad**28704**  
**EPA - SPCC - U - O - UL**

Código Aduanero

**94032090**

Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Capacidad de almacenamiento (L)	Cubetas de retención	Estanterías	Carga admisible	Peso (kg)
28704	229 x 2464 x 635		800		0	0	35.00


**TRIONYX**  
DEPUIS 1971

 138 avenue Henri BARBUSSE 92140 Clamart (France)  
 Tél. 01 46 45 80 00 - Fax : 01 46 45 85 98  
 www.trionyx.fr - trionyx[at]trionyx.fr