

1619R



Flaconi di sicurezza con Viton O-ring – 0,4 L

CONFORMIDAD

ADVANTAGES

- Retention equipment
- Storage of flammables and toxics

SINO TAMBIÉN...

- Ausencia de fugas cuando el frasco se vuelca
- Sealing cap with nitrile rubber joint
- Available in 2 colors (red or yellow)
- Available in different sizes and capacities
- Accurate dosing
- Double opening and sealing O-ring and spring
- Cierre automático
- Lid opening by finger pressure



CARACTERISTICAS TECNICAS

Modelo	1619R			Código Aduanero	73102910		
Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Capacidad de almacenamiento (L)	Cubetas de retención	Estanterías	Carga admisible	Peso (kg)
1619R	222 x 65 x 65		0		0	0	0.10