

1633R



## Flacons di sicurezza con Viton O-ring – 1 L

### CONFORMIDAD

### ADVANTAGES

- Retention equipment
- Storage of flammables and toxics

### SINO TAMBIÉN...

- Ausencia de fugas cuando el frasco se vuelca
- Sealing cap with nitrile rubber joint
- Available in 2 colors (red or yellow)
- Available in different sizes and capacities
- Accurate dosing
- Double opening and sealing O-ring and spring
- Cierre automático
- Lid opening by finger pressure



### CARACTERISTICAS TECNICAS

Modelo

1633R

Código Aduanero

73102910

Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Capacidad de almacenamiento (L)	Cubetas de retención	Estanterías	Carga admisible	Peso (kg)
1633R	280 x 80 x 80		1		0	0	0.20