

10608



## Recipienti di sicurezza in acciaio stagnato per prodotti infiammabili - Tipo Laboratorio - 9,5 L

### CONFORMIDAD

FM

### ADVANTAGES

- NFPA and OSHA compliant
- Conforme a la norma FM
- Storage of flammables and toxics

### SINO TAMBIÉN...

- FM, NFPA, OSHA and UL approved
- Available in different sizes and capacities
- Available in different colors



### CARACTERISTICAS TECNICAS

Modelo  
Conformidad**10608**  
**FM**

Código Aduanero

**73102900**

Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Capacidad de almacenamiento (L)	Cubetas de retención	Estanterías	Carga admisible	Peso (kg)
10608	292 x 298 x 298		10		0	0	4.00