

CORG203



NFX  
15-211




## Filtro ha carbone attivo polivalente - Compatible MF3-BE

### CONFORMIDAD

### ADVANTAGES

- Specific filtration of the chemicals
- NF X 15-211 compliant
- Storage of flammables, corrosives and toxics



### CARACTERISTICAS TECNICAS

Modelo	CORG203			Código Aduanero	38021000		
Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Capacidad de almacenamiento (L)	Cubetas de retención	Estanterías	Carga admisible	Peso (kg)
CORG203	100 x 565 x 405		0		0	0	13.00