

TROR



NFX
15-211



Zoppico di 10 tubi reagenti

CONFORMIDAD

ADVANTAGES

- Convenient for the measure of the filter's saturation
- NF X 15-211 compliant
- Storage of flammables, corrosives and toxics



CARACTERISTICAS TECNICAS

Modelo	TROR	Código Aduanero	32041600
Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Peso (kg)
TROR			1.00