

TROR



NFX  
15-211



## Zoppico di 10 tubi reagenti

### CONFORMIDAD

### ADVANTAGES

- Convenient for the measure of the filter's saturation
- NF X 15-211 compliant
- Storage of flammables, corrosives and toxics



### CARACTERISTICAS TECNICAS

Modelo	TROR	Código Aduanero	32041600
Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Peso (kg)
TROR			1.00