



Armadio di sicurezza bicolore porte blu - 300L

CONFORMITÀ

Articles R.5132-66 et R.5132-68 Code de la santé publique

ADVANTAGES

- Compliant to the Public Health Code
- Conforma alle raccomandazioni dell INRS
- Storage of flammables, corrosives and toxics

MA ANCHE...

- Compliant to the Public Health Code articles R. 5132-66 and R. 5132-68
- Conforma alle raccomandazioni dell INRS
- Allows the storage of your chemicals safely
- Affordable equipment to store your chemicals



CARATTERISTICHE TECNICHE

Modello	AZ300B	Chiusura porta	chiusura manuale
Conformità	Articles R.5132-66 et R.5132-68 Code de la santé publique	Dimensioni di spedizione HxWxD (mm)	2150 x 1200 x 800
Materiale	acciaio	Ventilazione	orifizio di ventilazione alto e bassi sulle porte
Colore porte	blu RAL 5010	Numero di compartimenti	1
Numero di porte	2	Vaschetta di contenimento	1
Tipo porte	porte battenti	Codice HS	94032080
Cerniere	cerniere a destra e sinistra		

COMPLEMENTARY ACCESSORIES - VENTILATION AND FILTRATION BOXES

- Opzione taglio copertura per raccordo (Réf : AZVENT)
- Cassettone di ventilazione in acciaio (Réf : CDV-A)
- Kit di raccordo per cassettoni (Réf : KRC)

OPTIONS

- Carabottino per armadio AZ110 ed AZ300 (Réf : CE300)
- Mensola/bac per armadio AZ110 ed AZ300 (Réf : ESE300)

Modello	Dimensioni esterne HxWxD (mm)	Dimensioni interne HxWxD (mm)	Capacità di stoccaggio (L)	Vaschetta di contenimento	Ripiani	Capacità di carico	Peso (kg)
AZ300B	1950 x 930 x 500	1840 x 925 x 453	300	1	3	80	90.00