

1619R



## Flaconi di sicurezza con Viton O-ring – 0,4 L

### CONFORMITÀ

### ADVANTAGES

- Self-closing doors
- Retention equipment
- Storage of flammables and toxics

### MA ANCHE...

- Assenza di fughe quando il flacone si rovescia
- Sealing cap with nitrile rubber joint
- Available in 2 colors (red or yellow)
- Available in different sizes and capacities
- Accurate dosing
- Double opening and sealing O-ring and spring
- Chiusura automatica
- Lid opening by finger pressure



### CARATTERISTICHE TECNICHE

Modello	1619R			Codice HS	73102910		
Modello	Dimensioni esterne HxWxD (mm)	Dimensioni interne HxWxD (mm)	Capacità di stoccaggio (L)	Vaschetta di contenimento	Ripiani	Capacità di carico	Peso (kg)
1619R	222 x 65 x 65		0		0	0	0.10