

L137

ISO
3864

EN
15154-1

EN
15154-2



Aérateur 22x100 avec tamis

CONFORMITÀ

ADVANTAGES

- ISO 3864 compliant
- First-aid in case of chemical projection



CARATTERISTICHE TECNICHE

Modello

L137

Codice HS

73242921

Modello	Dimensioni esterne HxWxD (mm)	Dimensioni interne HxWxD (mm)	Flusso di lavaggio oculare (L/min)	Flusso doccia (L/min)	Peso (kg)
L137			0	0	0.60